

Astor Motion G Tilt in Space Reclining Bath with shower gantry



Highlights:

- Free standing with height and reclining adjustment - simple electronic operation
- Extra large lip-free door opening and low entry height of 49cm: easy self entry
- Large footwell (98 ltr) with auto pre-fill
- Fastest drain and speedy fill time
- Interior moulded seat for extra comfort and support
- WRAS approved shower gantry

The Astor Motion G reclining bath is an iconic bath offering the very latest style combined with a revolutionary approach to bathing. Our specialist design team has worked very closely with professionals in the care sector to combine a contemporary appearance with practical performance.

Despite its advanced mechanism it takes just one carer to operate the Astor Motion G, easily adjusting the height and tilt of the bath to allow easy access for the bather. The side door gives a large lip-free aperture making access easier for wheelchair users, self-bathers or to be used in combination with a hoist. Space under the bath also allows hoists or wheelchairs to be positioned very close to the bath for a smoother transfer.

The large footwell (98 litres) can be filled before entry to allow a bather to have enough water to be immersed as soon as the bath is tilted back. The autofill mode is an additional aid to the carer to allow attention to be given elsewhere whilst the bath is filling. This bath also boasts the quickest drain time (75 litres a minute) and a very speedy fill time (66 litres per minute), ensuring that the bather is immersed quickly as the bath fills and is not left in draining lukewarm water as the bath empties.

The Motion G comes with a WRAS approved shower gantry system with TMV3 temperature control and backflow protection to fluid category 5 level.



Standard Features:

- Maximum Working Load (MWL) 210 kg/33 st
- Height range 810 mm - 1010 mm (200 mm)
- External dimensions 1908 mm (L) x 935 mm (W)
- Tilt angle 23°
- Height adjustable (mains electric)
- Emergency battery backup
- WRAS approved shower gantry system
- Floor and ceiling hoist access
- Chrome grab handles
- Large lightweight door (damped safety hinged)
- Door handle with built in safety lock
- 98 ltr foot well (auto pre-fill options)
- 3 Year warranty (*Terms and conditions apply*)

Optional Extras:

- Air spa system
- Integrated Bluetooth music/sound system
- Colour changing LED light system
- Range of service & warranty packages.
Please ask for further details.

Product Specification Sheet:

Astor Motion G



Model Name:	Astor Motion G
Product Description	Height adjustable tilt in space reclining bath with shower gantry
Product Codes	02.108.00.000x 4 Motion G with TMV3 WRAS gantry & shower + handles
Height Range	1908mm (L) x 935mm (W)
Dimensions	810mm - 1010mm (200mm) / 23° Tilt angle

Features:

Height Adjustable	Mains Electric	23° Tilt Reclining Tub	98 ltr Foot Well	210kg/33st MWL	WRAS Gantry System

Options:

Range of Taps	Spa System	Bluetooth	LED Light System	Service Packages	Many More Options

Application:

Care Home	Supported Living	Multi-User Environment	Hospital

Technical Drawings: All dimensions in mm:

Note: A minimum of 250mm clearance is required on all sides of the Motion bath

